Flexible Reamer System

| | | | |
|----|---------------------------------------|----|---------------------------------------|
| EN | FLEXIBLE REAMER SYSTEM 2 | JA | フレキシブルリーマーシステム 72 |
| IT | SISTEMA DI ALESATORI FLESSIBILI 7 | ZH | 可屈式扩髓系统 78 |
| FR | SYSTÈME D'ALESAGE FLEXIBLE 13 | CS | FLEXIBILNÍ FRÉZOVACÍ SYSTÉM 83 |
| DE | FLEXIBLES FRÄSERSYSTEM 19 | PL | SYSTEM GIĘTKIEGO POSZERZACZA 89 |
| ES | SISTEMA DE FRESAS FLEXIBLES 25 | SL | SISTEM FLEKSIBILNIH SVEDROV 95 |
| DA | FLEXIBLE REAMER SYSTEM 31 | KO | 굴곡식 확공기 장치 101 |
| FI | TAIPUISA AVARRUSPORAJÄRJESTELMÄ 37 | AR | جهاز موسع الثقوب المرن 107 |
| NO | FLEKSIBELT BROTSJSYSTEM 43 | BG | ГЪВКАВА СИСТЕМА ЗА РАЗШИРЯВАНЕ 114 |
| NL | FLEXIBEL RUIMERSYSTEEM 48 | LT | LANKSTI PLĖSTUVO SISTEMA 120 |
| PT | SISTEMA DE FRESAS FLEXÍVEIS 54 | MS | SISTEM RIMER FLEKSIBEL 126 |
| SV | FLEXIBELT BROTSCHSYSTEM 60 | RO | SISTEM ALEZOR FLEXIBIL 132 |
| EL | ΕΥΚΑΜΠΤΟ ΣΥΣΤΗΜΑ ΔΙΕΥΡΥΝΤΗΡΑ 66 | TR | FLEXIBLE REAMER SİSTEMİ 138 |

ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone +39 045 6719000 - Fax +39 045 6719380

0123 PQFRS G 03/22 (0425544)

Flexible Reamer

System

EN

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online

Important information - please read prior to use

See also instruction leaflet PQRMD for reusable medical devices

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Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

ONLY USE THE FLEXIBLE REAMER SYSTEM WITH AN OLIVE-TIPPED GUIDE WIRE

NON STERILE

DESCRIPTION

ORTHOFIX instrumentation consists of medical devices and their accessories used in surgical procedures. The range includes manual surgical instruments and surgical instruments driven by power sources. ORTHOFIX devices should only be used by qualified personnel trained in the use of surgical instruments and the relevant surgical procedures. The Orthofix Flexible Reamer System comprises a range of intramedullary reamer heads in two different systems:

1. Tibia & Femur:
 - a. Monobloc 8,0 and 8,5mm reamer heads which are permanently attached to a flexible drive shaft;
 - b. Modular cutting heads from 9,0mm to 17,0mm in 0,5mm steps. These reamer heads connect to a flexible drive shaft that in turn is connected to a manual or powered driver. The cutting heads are locked in place when the reamer and drive shaft are passed over a 980mm x 3.0mm olive tipped guide wire.
2. Humerus: a set of monobloc flexible reamers from 7,0mm to 10,0mm in 0,5mm steps. They must be used over a 780mm x 2.0mm olive tipped guide wire.

Both Tibia & Femur and Humerus Systems have the following characteristics materials:

- Reamer heads, couplings and end fittings - stainless steel
- Flexible drive shaft (modular)- nickel titanium alloy
- Monobloc flexible reamers - head and base fittings of stainless steel; shaft of nickel titanium alloy

COUPLING

All the drive shafts, monobloc or modular, have AO quick connect couplings as standard. A Zimmer adapter is available in the kits.

INTENDED USE

The flexible reamer system is intended to clear the medullary canal of bone marrow and debris and to effectively enlarge the medullary canal of long bones for the acceptance of an intramedullary nail.

PRECAUTIONS

- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation.
- When an ORTHOFIX Flexible Reamer is used, in the tibia and femur a 3mm diameter 980mm long olive-tipped guide wire must be inserted first, using an image intensifier to cross a fracture site, and to check that the position of the guide wire is correct when fully inserted. In the humerus a 2.0mm diameter 780mm long olive tipped guide wire must be used in the same way.
- Check the position of the guide wires frequently using an image intensifier (fluoroscopy) to prevent unintended guide wire advancement and penetration into the surrounding tissues. Use the fluoroscope to monitor the advancement of the reaming head. Take care not to withdraw the guide wire when the reamer is extracted, and whenever there is any doubt re-check the position with fluoroscopy before passing the next reamer head.
- Prior to any use, guide wires should be checked and should be discarded if found to be scratched or bent.
- When using guide wires with instruments and/or cannulated screw systems clean the cannulated instruments intra-operatively to prevent the accumulation of bone debris in the cannulation. Clear the cutting flutes of bone debris on a regular basis.
- The width of the isthmus is determined by pre-operative X-ray examination. The instrument with the smallest possible diameter is used for initial reaming of the intramedullary canal. This will usually be 7.0mm in the humerus, 8.0mm in the tibia and 9.0mm in the femur. However individual bones vary and if the chosen reamer will not pass easily a smaller size should be selected. The reamer size should be increased in 0.5mm steps until the desired size is reached.

WARNINGS

- Care should be taken not to cut through surgical gloves when handling any sharp-edged surgical instrument and to take into account the risk of infection if a cut appears.
- Surgical instruments of this type should be inspected before use; if defects or cracks are obvious then the product should not be used.
- Do not use a cutting instrument that has dull edges.
- The reaming cutter heads **MUST** be used in the correct order, at 0.5mm intervals. Do not use a reaming system if there is a cutting head missing.

- It is important to avoid exerting too much pressure during reaming, to avoid the risk of fat embolism. Reaming should be performed with a light touch and gentle forward movement. If this is difficult, remove the reamer head and clean it, and if necessary pass the previous cutting head. Check that the reaming heads are being used in the correct order.
- The instrument should not be forced when resistance is encountered, and should always be withdrawn when it is not advancing. If jammed, the reamer can be freed by reversing the reamer driver direction and withdrawing from the canal. Excessive torque and/or stalling of the reamer may cause shaft fracture and fragmentation to occur.
- Penetration of the cortex has been reported with intramedullary reamers. An olive tipped guide wire should be used in all cases and fluoroscopy should be used to monitor the depth of cut. The guide wire is essential to lock the cutting head to the drive shaft, and the olive is used to help extract a jammed cutting head.
- To prevent over-reaming of the anterior tibial cortex at the entry portal, ensure that the reamer head is fully inserted into the canal before applying power.
- Instrument failure or a reduced useful life is more likely to occur with misuse. Small scratches can have a significant effect on the fatigue life or strength of the product.
- Even with correct use, care and maintenance surgical instruments should not be expected to last indefinitely. This is especially true for cutting instruments (e.g., drills, reamers, taps and rasps) and driving instruments (e.g., mallets, impactors, pins and extractors). Such instruments are subjected to high loads and/or impact forces. Under such conditions breakage can occur, particularly when the instrument is corroded, damaged or scratched.

STERILE & NON-STERILE PRODUCTS

Orthofix devices are provided as STERILE or NON-STERILE and they are labeled as such.

In the case of STERILE products, product integrity, sterility and performance are assured only if the packaging is undamaged. Do not use if packaging is compromised, unintentionally opened or if a component is believed to be faulty, suspect or damaged.

The products supplied NON-STERILE require cleaning, disinfection and sterilization prior to use according to procedures reported in the following instructions.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use.

DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.

Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.
5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.

8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

CLEANING

General considerations

Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5 % anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300 Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution to displace trapped air; ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three times with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using absorbent not shedding cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.
6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
8. Arrange medical devices to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 minutes;

- b. Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water for 10 minutes at 55°C;
- c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 minutes;
- d. Final rinsing with deionized water for 3 minutes;
- e. Thermal disinfection at least 90°C or 194° F (max 95 °C or 203° F) for 5 minutes or until A0=3000 is reached. The water used for thermal disinfection must be purified.
- f. Drying at 110 °C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.

The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.

10. Select and start a cycle according to the recommendations of the washer manufacturer.
11. On completion of the cycle, ensure that all stages and parameters have been achieved.
12. Wearing protective equipment unload the washer disinfector when it completes the cycle.
13. If necessary, drain excessive water and dry by using absorbent not shedding cloth.
14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use.

All functional checks and inspections described below also cover the interfaces with other instruments or components.

The failure modes below may be caused by end of life of the product, improper use or improper maintenance.

Orthofix does not typically specify the maximum number of uses for re-usable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use are the best methods of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as per the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED. In particular, if cracks or scratches are visible on the titanium-nickel alloy (Nitinol) shaft, the instrument must not be used under any circumstances and must be absolutely replaced by a new one. Tests prove that these defects are the cause of shaft breakage, which may take place during use in the patient's medullary cavity with the risk of not being able to remove all the fragments.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use.

Specific instructions may be available for some product codes. These instructions are linked to the product code and are available on a dedicated Orthofix website.

Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS).The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory. In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray.

Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products.

Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization.

Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

| Steam sterilizer type | Gravity | Pre-vacuum | Pre-vacuum | Pre-vacuum |
|------------------------------|-------------------|-----------------|-------------------|-----------------|
| Notes | Not for use in EU | - | Not for use in US | WHO guidelines |
| Minimum exposure Temperature | 132° C (270° F) | 132° C (270° F) | 134° C (273° F) | 134° C (273° F) |
| Minimum exposure Time | 15 minutes | 4 minutes | 3 minutes | 18 minutes |
| Drying Time | 30 minutes | 30 minutes | 30 minutes | 30 minutes |
| Number of pulses | N/A | 4 | 4 | 4 |

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0.5%

IMPORTANT STATEMENT

- It is strictly prohibited to carry out any modification whatsoever on an ORTHOFIX instrument. If this recommendation is not followed, ORTHOFIX disclaims any liability for any subsequent consequential event.
- The product is marketed for the specific indications described in the labelling. The use of this product for other than its intended purpose is either contra-indicated or is without evidence to support the safety and effectiveness of such use. For the information of individuals and institutions contemplating the use of this product for other than labelled indications (i.e., off-label use), such use may be experimental and may be the subject of restrictions under applicable laws and regulations.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

| Symbol | Description | |
|---|--|--|
|   | Consult instructions for use or consult electronic instructions for use | CAUTION: Consult instructions for use for important cautionary information |
|  | Single Use. Do not re-use | Orthofix note: discard appropriately after use (treatment) on patient |
|  | Sterile. Sterilized using irradiation | |
|  | NON STERILE | |
|   | Catalogue number | Batch code |
|  | Use-by date (year-month-day) | |
|  | CE marking in conformity to applicable European Medical Device Directives/Regulations | |
|   | Date of manufacture | Manufacturer |
|  | Do not use if package is damaged and consult instructions for use | |
| Rx Only | Caution: Federal law (USA) restricts this device to sale by or on the order of a physician | |

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Orthofix Srl