
ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone 0039-0456719000 - Fax 0039-0456719380

**EN****Instruction leaflet - Please read prior to use**

TrueLok

External fixation system TrueLok Screws

Consult instructions for Use

RX Only

Important information for the surgeon performing the operation

Read instructions before use

INFORMATION TO USE TrueLok SCREWS (HALF PINS)

TrueLok HALF PINS SHOULD BE USED ONLY WITH THE TrueLok SYSTEM FOR SPECIFIC INFORMATION ON THE TrueLok SYSTEM REFER TO PQ TLK.

It is essential that special care be taken to ensure that the screws do not penetrate joints or growth cartilage in pediatric patients. The length of both the entire screw and the threaded section must be selected in accordance with the dimensions of the bone and soft tissues.

Screws should always and only be inserted using a driver or T-wrench, regardless of whether or not pre-drilling has been performed. It is important to apply moderate force to penetrate the first cortex. Diaphyseal bone screws should always be inserted in the centre of the bone axis, in order to avoid weakening. In all cases, the surgeon should be mindful of the amount of torque required to insert the screw. Should this appear more difficult than usual, it is recommended the screw be removed and cleaned and the bone drilled once more using the drill bit, even if this has already been used.

Excessive penetration of the second cortex must be avoided because this may cause soft tissue damage. Bone screws must never be inserted so far that the smooth shaft penetrates the first cortex, in order to avoid the risk of bone damage.

The diameter of the screw must be selected in accordance with the diameter of the bone: for a bone diameter greater than 20 mm, 6.0 mm, 6.0 - 5.0 or 6.0 - 5.6 mm should be used; for a bone diameter of between 15 and 20 mm, 5.0 mm screws should be used and for a bone diameter between 10 mm and 15 mm, 4.0 mm bone screws are recommended.

Self-drilling screws with a thread diameter of 5.0 mm or more must never be inserted with a power tool, and should be inserted manually or with the use of a hand drill. Self-drilling screws with a smaller thread diameter may be inserted either manually or with the use of a low-speed electric drill.

TrueLok screws are not tapered and have a thread that has the same diameter throughout.

TrueLok screws are available in 4, 5, and 6mm diameters, each having a total length of 180mm and a thread length of 55mm. If the screw has to be inserted in areas wider than 50mm, XCaliber bone screws should be used. TrueLok screws are self-drilling, self-tapping, and have a standard quick connect end for ease of use in the operating room.

The following steps outline the recommended method of screw insertion and fixation:

1. A half pin fixation bolt is inserted into an appropriate hole of the ring.
The half pin fixation bolt acts as a guide for half pin insertion
2. A K-wire is passed through the hole of the fixation bolt along the intended direction of the half pin insertion. Make a stab incision in the skin at this level
3. Use scissors or hemostats to make a track through the soft tissues to the bone with blunt dissection only.
4. The desired half pin is inserted through the fixation bolt and soft tissue track, and drilled into the first bony cortex. Insertion through the second cortex should be done manually using the dedicated Orthofix T-wrench.
5. The half pin fixation bolt is then secured firmly to the external support with a nut.

If TrueLok screws are cut to length, screw covers should be used to protect the cut end. TrueLok screws which have been cut to length could be removed manually with the Universal Chuck with

T-Handle (17955) or with a power drill.

Care should be taken of the soft tissues during screw or wire insertion.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

SYMBOL

Consult instructions for Use

Single use. Do not reuse

Sterilised with ethylene oxide

Sterilised by irradiation

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Catalogue number / Lot number

Expiry date (year-month)

CE marking in conformity to MDD 93/42/ECC as amended by 2007/47/EC

Date of manufacture / Manufacturer

Do not use if package is opened or damaged

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