Contours VPS™

Volar Plating System

EN	CONTOURS VPS™ VOLAR PLATING SYSTEM 2	JA	CONTOURS VPS 掌側プレートシステム 74
IT	CONTOURS VPS SISTEMA DI PLACCHE VOLARI 8	ZH	CONTOURS VPS 掌侧钢板系统 80
FR	CONTOURS VPS PLAQUE PALMAIRE 14	CZ	CONTOURS VPS SYSTÉM VOLÁRNÍHO DLAHOVÁNÍ 86
DE	CONTOURS VPS VOLARES PLATTENSYSTEM 20	PL	CONTOURS VPS SYSTEM PŁYTEK DŁONIOWYCH 92
ES	CONTOURS VPS SISTEMA DE PLACA VOLAR 26	SL	CONTOURS VPS SISTEM VOLARNE PLOŠČE 98
DA	CONTOURS VPS VOLÆRT PLADESYSTEM 32	KO	CONTOURS VPS VOLAR 플레이팅 시스템 104
FI	CONTOURS VPS VOLAR-LEVYJÄRJESTELMÄ 38	AR	CONTOURS VPS جهاز التصفيح الراحي 110
NO	CONTOURS VPS VOLARPLATESYSTEM 44	BG	CONTOURS VPS СИСТЕМА ЗА ВОЛАРНИ ЗАКЛЮЧВАЩИ ПЛАСТИНИ 117
NL	CONTOURS VPS VOLAIR PLAATSYSTEM 50	LT	CONTOURS VPS DELNINIO FIKSAVIMO PLOKŠTELĖMIS SISTEMA 123
PT	CONTOURS VPS SISTEMA DE PLACA VOLAR 56	MS	CONTOURS VPS ESISTEM PEMPLATAN VOLAR 129
SV	CONTOURS VPS SYSTEM FÖR VOLARPLATTA 62	RO	CONTOURS VPS SISTEM DE FIXARE VOLARĂ CU PLACĂ 135
EL	CONTOURS VPS ΣΥΣΤΗΜΑ ΠΑΛΑΜΙΑΙΑΣ ΣΤΕΡΕΩΣΗΣ 68	TR	CONTOURS VPS VOLAR PLAKA SİSTEMLİ 141

ORTHOFIX®
ORTHOFIX Srl
Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
Telephone +39 045 6719000 - Fax +39 045 6719380

PQVPS H 02/22 (0424728)

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online

Important information - please read prior to use

See also instruction leaflet PQRMD for reusable medical devices

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Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

INDICATIONS FOR USE

The Orthofix Contours VPS, Volar Plating System, is intended for the volar fixation of fractures and osteotomies involving the distal radius.

DESCRIPTION

The Contours VPS is an anatomically contoured, Delta-shaped plate intended for volar applications to the distal radius. The benefits to utilizing a volar approach for fixation of distal radius fractures and osteotomies are that the dorsal tissues are preserved and it may reduce the chances of tendon ruptures occurring.

The Orthofix Contours VPS Plate features an angled head with two rows of holes for placement of screws distally on the radius to secure bone fragments and provide stabilization. If a bone graft is required, access to the site is made possible by the triangular opening in the plate that spans from the plate head to the distal portion of the shaft. The shaft of the Plate is placed proximally on the distal radius and offers a variety of K-wire and screw placement options for secure fixation to the external surface of the bone. The device provides the stabilization and fixation necessary in the treatment of distal radius fractures and osteotomies to allow early functional use of the hand.

CONTRAINDICATIONS

- Cases in which there is active infection local to the operative site.
- Any patient having inadequate tissue coverage over the operative site or where there is inadequate bone stock, bone quality, or anatomical definition.
- Patients with metabolic bone disorders such as severe osteopenia, osteoporosis, osteomyelitis, poorly controlled diabetes mellitus and Paget's disease.
- Cases in which the Contours VPS Plate would cross joint spaces or open epiphyseal growth plates such as in skeletally immature patients.
- Any time implant utilization would interfere with anatomical structures or physiological performance.
- Suspected or documented metal allergy or intolerance.
- Any patient unwilling or incapable of following post-operative care instructions.

WARNINGS AND PRECAUTIONS

- Bone plates and screws must not be reused. Reuse could result in injury or reoperation due to breakage or infection. Do not attempt to resterilize implants that come in contact with body fluids.
- The Contours VPS Plate has not been evaluated for safety and compatibility in the MR environment. Contours VPS Plate has not been tested for heating or migration in the MR environment.
- Smoking, chronic steroid use and the use of other anti-inflammatory drugs have been determined to affect bone healing and could potentially have an adverse effect on the bone regenerate during fracture healing.
- The Contours VPS Plate is available in a RIGHT and LEFT version, which are anatomically shaped. It is therefore important to select the correct plate. DO NOT use a plate labelled "RIGHT" for a left wrist and vice versa.
- It is essential that proper operative technique be followed for implantation. Please refer to the appropriate Operative Technique Guide.
- Care should be used in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Scratches may reduce the functional strength of the component.
- Examine all Contours VPS components carefully prior to use to assure there is no damage to any of the components. If you believe a component to be faulty, damaged, or suspect **DO NOT USE**.
- The Contours VPS Plate and Bone Screws should not be combined with components from another manufacturer or system.
- Different metal types should not be used together. Combining implants of dissimilar metals may result in galvanic corrosion.
- The Steri-Tray is for sterilization of Contours VPS components **ONLY**.

POSSIBLE ADVERSE EFFECTS

- Bending or breaking of the plate or loosening, bending, or breaking of the bone screws in situ.
- A secondary surgical procedure may be necessary to replace or remove the bone plate and screws.
- Inadequate tissue surrounding the implant may cause localized irritation and/or pain.
- Scar formation possibly causing pain and/or neurological compromise around nerves.
- Non-union, delayed union, or mal-union could occur.
- Bone loss or decrease in bone density, possibly caused by stress shielding.
- Intractable pain.

- Complications associated with metal sensitivity.
- Intrinsic risks associated with anesthesia and surgery. Hemorrhage, hematoma, seroma, embolism, edema, stroke, excessive bleeding, phlebitis, wound necrosis, wound infection, or damage to blood vessels.

IMPORTANT

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons, or device failures that require surgical re-intervention to remove or replace the plate or bone screws. Proper implantation and patient selection will greatly affect the results.

MATERIALS

The Contours VPS Plate and Bone Screws are currently available in titanium alloy. Please refer to the individual product labels for component materials.

It is important that only components of similar materials are used together to eliminate the risk of galvanic corrosion.

STERILE AND NON-STERILE PRODUCT

The Volar plating system manufactured by Orthofix Srl is Gamma-sterile.

- A set of templates has been developed specifically for the sterile version of the Volar plating system. The templates are provided in non-sterile configuration.
- Please review the product label to determine the sterility of each device.
- Contents of package are STERILE unless package is opened or damaged. Do not use if package is opened or damaged.
- Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures. Please refer to section below.
- Product integrity and performance are assured only if packaging is undamaged.

INSTRUCTIONS FOR PROCESSING AND REPROCESSING

These reprocessing instructions have been written in accordance with ISO17664 and have been validated by Orthofix in compliance with international standards. It is the responsibility of the healthcare facility to ensure that the reprocessing is performed in accordance with instructions provided.

Warnings

- Devices labeled "FOR SINGLE USE ONLY" can be reprocessed multiple times before their first clinical use but must not be reprocessed for re-use in a clinical setting.
- Single use devices MUST NOT BE REUSED, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications. Please refer to the device label to identify single or multiple use and/or cleaning and re-sterilization requirements.
- The personnel that works with contaminated medical devices must follow safety precautions as per the procedure of the healthcare facility.
- Cleaning and disinfection solutions with a pH 7-10.5 are recommended. Cleaning and disinfection solutions with a higher pH should be used according to the material compatibility requirements stated on the detergent technical data sheet.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions MUST NOT be used.
- Contact with saline solutions should be minimized.
- Complex devices such as those with hinges, lumens or mated surfaces must be thoroughly manually pre-cleaned before automated washing to remove soiling that accumulates in recesses. If a device needs particular care in pre-cleaning, a product specific IFU is available on the Orthofix website, which is accessible using the data matrix reported on the product labelling.
- DO NOT use metal brushes or steel wool.

Limitations on reprocessing

- Repeated reprocessing has minimal effect on devices for which reprocessing is allowed.
- End of life is normally determined by wear and damage due to use.
- Products labeled for Single Use only MUST NOT be reused regardless of any reprocessing in a clinical setting.

POINT OF USE

Reprocess the reusable medical devices as soon as is reasonably practicable to minimize the drying of soil and residuals. For optimal results, instruments should be cleaned within 30 minutes of use. DO NOT use a fixating detergent or hot water because these can cause the fixation of residue.

CONTAINMENT AND TRANSPORTATION

Cover contaminated instruments during transportation to minimize the risk of cross contamination. All used surgical instruments must be regarded as contaminated. Follow the hospital protocols for handling contaminated and bio-hazardous materials. Handling, collection and transportation of used instruments must be strictly controlled to minimize any possible risks to patient, personnel and any area of the healthcare facility.

PREPARATION FOR CLEANING

This procedure can be omitted in case of direct subsequent manual cleaning and disinfection. In the case of a highly contaminated reusable medical device, before starting an automatic cleaning process, a pre-cleaning and a manual cleaning (described below) are recommended.

Manual Pre-cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient detergent solution. Orthofix recommends the use of a slightly alkaline enzymatic detergent solution based on a detergent containing <5% anionic surfactants and enzymes prepared using deionized water.
4. Carefully immerse the component in the solution to displace trapped air.

5. Scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Clean the single components using an ultrasonic device in a degassed cleaning solution.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed. Use a syringe in case of lumens or cannulations.
11. Remove item from rinse water and drain.
12. Carefully hand-dry using absorbent not shedding cloth.

CLEANING

General considerations

Orthofix provides two methods of cleaning: a manual method and an automated method. Wherever applicable, the cleaning phase should start immediately after the pre-cleaning phase to avoid soil drying. The automated cleaning process is more reproducible and therefore more reliable, and the staff is less exposed to the contaminated devices and to the cleaning agents used. Staff shall comply with the safety precautions and procedures of the healthcare facility regarding the use of protective equipment. In particular, staff should take note of the instructions provided by the cleaning agent manufacturer for correct handling and use of the product. Observe all instructions provided by the detergent manufacturer regarding the immersion time of the device in the cleaning agent/disinfectant and its concentration. The quality of the water used for diluting cleaning agents and for rinsing medical devices should be carefully considered.

Manual cleaning

1. Wear protective equipment in compliance with the safety precautions and procedures of the healthcare facility.
2. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
3. Fill the receptacle with sufficient cleaning solution. Orthofix recommends the use of a slightly alkaline enzymatic cleaning solution.
4. Carefully immerse the component in the solution to displace trapped air; ensure that the cleaning solution reached all surfaces, including holes or cannulations.
5. Thoroughly scrub the device in the cleaning solution with a soft bristle nylon brush until all visible soiling is removed. Use a soft bristle nylon brush to remove residuals from lumens, with a twisting motion on rough or complex surfaces.
6. Rinse cannulations at least three times with cleaning solution using a syringe. Never use metal brushes or steel wool.
7. Remove the device from the cleaning solution.
8. Brush the single components in running tap water.
9. Put single components in an ultrasonic device with degassed cleaning solution at 2% for 10 minutes. Orthofix recommends the use of a detergent solution based on a detergent containing < 5 % anionic surfactants, non-ionic surfactants and enzymes, prepared using deionized water. Orthofix recommends using an ultrasound frequency of 35kHz, power = 300 Weff, time 15 minutes. The use of other solutions and parameters shall be validated by the user and the concentration shall be in compliance with the detergent manufacturer's technical datasheet.
10. Rinse the components in purified sterile water until all traces of cleaning solution are removed.
11. Rinse the cannulations, rough or complex surfaces at least three times with purified sterile water. When cannulations are present use a syringe to facilitate this step.
12. Remove item from rinse water and drain.
13. If, after completion of the cleaning steps, some encrusted soil remained on the device, the cleaning steps must be repeated as described above.
14. Carefully hand-dry using absorbent not shedding cloth.

Manual disinfection

1. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
2. Fill the receptacle with sufficient disinfectant solution. Orthofix recommends the use of a 6% hydrogen peroxide solution for 30 minutes prepared using water for injection.
3. Carefully immerse the component in the solution to displace trapped air; ensure that the disinfectant solution reached all surfaces, including holes or cannulations.
4. Rinse cannulations, rough or complex surfaces at least three times with disinfectant solution. Use a syringe filled with disinfectant solution to rinse cannulations.
5. Remove the items from the solution and drain.
6. Soak in water for injection (WFI) to remove traces of disinfectant solution.
7. Rinse the cannulations at least three times with a syringe (filled with WFI).
8. Remove item from rinse water and drain.
9. Repeat the rinsing procedure as described above.
10. Carefully hand-dry using absorbent not shedding cloth.
11. Visually inspect and repeat manual cleaning and disinfection if necessary.

Automatic cleaning and disinfection using washer-disinfector

1. Perform a pre-cleaning if necessary due to the contamination of the device. Take special care when the items to be cleaned contain or have:
 - a. Cannulations
 - b. Long blind holes
 - c. Mating surfaces
 - d. Threaded components
 - e. Rough surfaces
2. Use a washer-disinfector in compliance with EN ISO 15883 that is properly installed, qualified and regularly subjected to maintenance and testing.
3. Ensure that the cleaning receptacle is clean and dry; no visible foreign material can be present.
4. Ensure that the washer-disinfector and all services are operational.
5. Load the medical devices into the washer-disinfector. Place heavier devices in the bottom of the baskets. Products must be disassembled before placing them in the baskets according to the specific instructions provided by Orthofix. Wherever possible, all parts of disassembled devices should be kept together in one container.

6. Connect cannulations to the rinsing ports of the washer-disinfector. If no direct connection is possible, locate the cannulations directly on injector jets or in injector sleeves of the injector basket. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
 7. Avoid contact between devices because movement during washing could cause damage to devices and the washing action could be compromised.
 8. Arrange medical devices to locate the cannulations in a vertical position and so blind holes incline downwards to promote the leakage of any material.
 9. Use approved thermal disinfection program. When using alkaline solutions, a neutralizer must be added. Orthofix recommends that cycle steps are at least as follows:
 - a. Pre-cleaning for 4 minutes;
 - b. Cleaning with the appropriate solution. Orthofix recommends the use of a detergent solution based on a detergent containing <5% anionic surfactants, non-ionic surfactant and enzymes, prepared using deionized water for 10 minutes at 55°C;
 - c. Neutralization with basic neutralizing agent solution. Orthofix recommends the use of a detergent solution based on citric acid, concentration 0.1% for 6 minutes;
 - d. Final rinsing with deionized water for 3 minutes;
 - e. Thermal disinfection at least 90°C or 194°F (max 95 °C or 203°F) for 5 minutes or until A0=3000 is reached. The water used for thermal disinfection must be purified.
 - f. Drying at 110 °C for 40 minutes. When the instrument has a cannulation, an injector should be used to dry the internal part.
- The suitability of other solutions, concentration, time and temperature shall be checked and validated by the user following the detergent manufacturer's technical datasheet.
10. Select and start a cycle according to the recommendations of the washer manufacturer.
 11. On completion of the cycle, ensure that all stages and parameters have been achieved.
 12. Wearing protective equipment unload the washer disinfecter when it completes the cycle.
 13. If necessary, drain excessive water and dry by using absorbent not shedding cloth.
 14. Visually inspect each device for remaining soil and dryness. If soil remains repeat the cleaning process as described above.

MAINTENANCE, INSPECTION AND FUNCTION TESTING

The following guidelines shall be applied to all Orthofix instruments that are labeled for multiple use.

All functional checks and inspections described below also cover the interfaces with other instruments or components.

The failure modes below may be caused by end of life of the product, improper use or improper maintenance.

Orthofix does not typically specify the maximum number of uses for re-usable medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional testing of the device before use are the best methods of determining the end of the serviceable life for the medical device. For sterile devices the end of life has been defined, verified and specified with an expiration date.

The following general instructions apply to all Orthofix products:

- All instruments and product components must be visually inspected under good light for cleanliness. If some areas are not clearly visible, use a 3% hydrogen peroxide solution to detect the presence of organic residuals. If blood is present, bubbling will be observed. After the inspection, the device shall be rinsed and drained as per the instruction given above.
- If visual inspection evidences that the device was not properly cleaned, repeat the cleaning and disinfection steps or discard the device.
- All instruments and product components must be visually inspected for any signs of deterioration that may cause failure during use (such as cracks or damage to surfaces) and functions tested before being sterilized. If a component or instrument is believed to be faulty, damaged or suspect, it must NOT BE USED.
- Products that show excessive fading of marked product code, UDI and lot, thus preventing clear identification and traceability, must NOT BE USED.
- Cutting instruments must be checked for sharpness.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate hinges and moving parts with an oil that does not interfere with steam sterilization as per manufacturer's instructions before sterilization. Do not use silicone based lubricant or mineral oil. Orthofix recommends the use of a highly purified white oil composed by paraffinum liquidum of food and pharmaceutical grade.

As a general preventive action Orthofix recommends following the instructions in the operative technique to avoid damages related to incorrect use.

Specific instructions may be available for some product codes. These instructions are linked to the product code and are available on a dedicated Orthofix website.

Moreover, it is important to follow the cleaning procedure suggested by Orthofix to avoid damages related to incorrect handling.

PACKAGING

In order to prevent contamination after sterilization Orthofix recommends using one of the following packaging systems:

- a. Wrap in compliance with EN ISO 11607, suitable for steam sterilization, and appropriate to protect the instruments or trays contained from mechanical damage. Orthofix recommends the use of a double wrap constituted of tri laminate non woven fabrics made of spunbond polypropylene and melt blown polypropylene (SMS).The wrap shall be resistant enough to contain devices up to 10kg. In the USA, a FDA cleared sterilization wrap must be used and compliance with ANSI/AAMI ST79 is mandatory In Europe, a sterilization wrap in compliance with EN 868-2 may be used. Fold the wrap to create a sterile barrier system according to a process validated as per ISO 11607-2.
- b. Rigid sterilization containers (such as Aesculap JK series rigid sterilization containers). In Europe, a container in compliance with EN 868-8 may be used. Do not include additional systems or instruments in the same sterilization container.

Every other sterile barrier packaging not validated by Orthofix must be validated by the individual healthcare facility according to instructions from the manufacturer. When equipment and processes differ from those validated by Orthofix the healthcare facility should verify that sterility can be achieved using parameters validated by Orthofix. Do not include additional systems or instruments in the sterilization tray. Note that sterility cannot be guaranteed if the sterilization tray is overloaded. The total weight of a wrapped instrument tray should not exceed 10kg.

STERILIZATION

Steam sterilization according to EN ISO 17665 and ANSI/AMMI ST79 is recommended. Gas plasma, dry heat and EtO sterilization MUST BE avoided as they are not validated for Orthofix products.

Use a validated, properly maintained and calibrated steam sterilizer. The steam quality must be appropriate for the process to be effective. Do not exceed 140°C (284°F). Do not stack trays during sterilization. Sterilize by steam autoclaving, utilizing a fractionated pre-vacuum cycle or gravity cycle according to the table below:

Steam sterilizer type	Gravity	Pre-vacuum	Pre-vacuum	Pre-vacuum
Notes	Not for use in EU	-	Not for use in US	WHO guidelines
Minimum exposure Temperature	132° C (270° F)	132° C (270° F)	134° C (273° F)	134° C (273° F)
Minimum exposure Time	15 minutes	4 minutes	3 minutes	18 minutes
Drying Time	30 minutes	30 minutes	30 minutes	30 minutes
Number of pulses	N/A	4	4	4

Orthofix recommends always using a pre-vacuum cycle for steam sterilization. The Gravity cycle was validated for wraps only, and it is suggested only when no other options are available. The Gravity cycle was not validated for sterilization in rigid containers.

STORAGE

Store the sterilized instrument in the sterilization packaging in a dry and clean environment at room temperature.

DISCLAIMER

The instructions provided above have been validated by Orthofix srl as being a true description for (1) processing a single-use device and a multiple use device for its first clinical use and (2) processing a multiple use device for its re-use. It remains the responsibility of the reprocessing officer to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes must be adequately recorded. Any deviation by the reprocessing officer from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences, and must also be appropriately recorded.

CLEANING AGENT INFORMATION

Orthofix used the following cleaning agents during validation of these processing recommendations.

These cleaning agents are not listed in preference to other available cleaning agents that may perform satisfactorily:

- For manual precleaning: Neodisher Medizym concentration 2%
- For manual cleaning: Neodisher Mediclean concentration 2%
- For automated cleaning: Neodisher Mediclean concentration 0.5%

SERVICE

For any information please contact your local distributor.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

RISKS DUE TO THE RE-USE OF "SINGLE USE" DEVICE

IMPLANTABLE DEVICE*

The "SINGLE USE" implantable device* of Orthofix is identified through symbol  reported on the product label.

After the removal from the patient, the implantable device* has to be discarded.

The re-use of implantable device* introduces contamination risks for users and patients.

The re-use of implantable device* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(*): Implantable device

Any device intended to be totally/partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is considered an implantable device.

NON IMPLANTABLE DEVICE

The "SINGLE USE" non implantable device of Orthofix is identified through symbol  reported on the label or indicated in the "Instructions For Use" supplied with the products. The re-use of "SINGLE USE" non implantable device cannot guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

Symbol	Description	
 	Consult instructions for use or consult electronic instructions for use	CAUTION: Consult instructions for use for important cautionary information
	Single use. Do not re-use	Orthofix note: discard appropriately after use (treatment) on patient.
STERILE 		STERILE. Sterilised by irradiation
		NON STERILE
REF 	LOT 	Catalogue number
		Use-by date (year-month-day)
		CE marking in conformity to applicable European Medical Device Directives/Regulations
 	Date of manufacture	Manufacturer
		Do not use if package is damaged and consult instructions for use
Rx Only		Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
T		Titanium and its Alloys

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