

**Quick User's Guide**  
Software Version 1.2

<b>WORKFLOW</b>	<b>1</b>
<b>1. REQUEST A NEW ACCOUNT</b>	<b>3</b>
<b>2. LOGIN PAGE</b>	<b>5</b>
<b>3. HOME PAGE</b>	<b>7</b>
<b>4. PATIENT → ADD NEW PATIENT</b>	<b>7</b>
<b>5. PATIENT → LIST OF PATIENTS</b>	<b>9</b>
<b>6. CASES → LIST OF CASES</b>	<b>9</b>
<b>7. CASES → ADD NEW CASE</b>	<b>11</b>
<b>8. DEFORMITY PARAMATERS</b>	<b>13</b>
<b>9. FRAME PARAMETERS</b>	<b>15</b>
9.1 FRAME PARAMETERS - PRINT PAGE	17
<b>10. POSTOPERATIVE ASSESSMENT</b>	<b>19</b>
<b>11. END OF CORRECTION</b>	<b>21</b>
<b>12. SCHEDULE</b>	<b>23</b>
<b>13. PRESCRIPTION</b>	<b>25</b>
<b>14. REPORT</b>	<b>27</b>
<b>15. CHECKUP</b>	<b>29</b>
<b>16. CHANGE PASSWORD</b>	<b>31</b>

Orthofix wishes to thank the following surgeons for their contribution to the development of this User's Guide:


Franz Birkholtz, M.D.

Alexander Cherkashin, M.D.

Mikhail Samchukov, M.D.

William Terrell, M.D.

The following provides a sequential overview of the process for a case management with TL-HEX software



The screenshot shows the TL-HEX website homepage. At the top left is the TL-HEX logo with the text 'TRUELOK HEXAPOD SYSTEM'. At the top right is the ORTHOFIX logo. Below the logos is a large image of a hexapod system. A navigation bar below the image contains four links: 'PRODUCT INFORMATION', 'TRAINING', 'FAQ', and 'CONTACT US'. The main content area features the text 'The future of the hexapod system built on solid circular fixation experience'. To the right of this text are two sections: 'MEMBERS AREA' with a 'LOGIN' button (callout 1) and 'NEW USERS' with a 'Request an Account' button (callout 2). At the bottom of the page is a footer with the links 'PRIVACY POLICY | COOKIES | CONTACT US'.

**TL-HEX**  
TRUELOK HEXAPOD SYSTEM

**ORTHOFIX**

PRODUCT INFORMATION TRAINING FAQ CONTACT US

The **future** of the hexapod system built on solid circular fixation **experience**

MEMBERS AREA  
LOGIN 1

NEW USERS  
Request an Account 2

PRIVACY POLICY | COOKIES | CONTACT US


---

**1. Request a New Account**

---


**Description**

- |          |   |
|----------|---|
| <b>1</b> | Click [Login] to access the software if the account is already activated and the user has his/her own username and password |
| <b>2</b> | Click [Request an Account] button and follow the onscreen instructions to obtain a valid access to the TL-HEX software.     |
-



## 2. Login Page

	Description
<b>HEADER</b> - This section is reported in all the application pages	
1	Click the TL-HEX logo to go to the product website <a href="http://www.tlhex.com">www.tlhex.com</a>
2	Click the Orthofix logo to go to the Company website <a href="http://intl.orthofix.com">intl.orthofix.com</a>
3	[ Log In] - if the user is <u>not</u> logged in By clicking [Log In], the user is forwarded to the Login page.  [ Log Out] - if the user is logged in By clicking [Log out], the user will be logged out by the application and return to the homepage
<b>FOOTER</b> - This section is reported in all the application pages	
4	Click [Privacy Policy] to display the Privacy Policy statement
5	Click [EULA Policy] to review the End User Licence Agreement subscribed for the application
6	Click [Cookies] to review the Cookies statement
7	Click [Contact us] to find how to contact the TL-HEX Customer Care support
8	Click [Instruction for Use] to access the TL-HEX document area
<b>Login</b>	
9	Enter a valid username and password in the appropriate text fields <i>Please note that username and password are case sensitive.</i>
10	Click [Sign In] button to access to the TL-HEX software
11	Click the link <i>Forgot your Username or Password?</i> and follow the displayed instructions to recover username and/or password
12	Current Software version






**3. Home Page****MENU STRUCTURE** - Reported in all the application pages when user is logged in

	Description
1	Click [Home Page] to return to the application Homepage
2	Click [Cases] to access to the Cases Menu functionalities <ul style="list-style-type: none"> <li>• Select [List of Cases] to view all the cases (see screen N°.6 for detail)</li> <li>• Select [Add New Case] to start a new case (see screen N°.7 for detail)</li> </ul>
3	Click [Patients] to access to the Patients Menu functionalities <ul style="list-style-type: none"> <li>• Select [List of Patients] to view all the registered patients (see screen N°.5 for detail)</li> <li>• Select [Add New Patient] to start a new patient registration (see screen N°.4 for detail)</li> </ul> All cases are related to a patient, therefore a new patient who doesn't exist must be created prior to begin a new case.
4	Click [Account] to access to the Account Menu functionalities <ul style="list-style-type: none"> <li>• Select [Change password] to invoke the change password procedure (see screen N°.16 for detail)</li> </ul>

**4. Patient → Add New Patient**

	Description
1	Insert the Patient ID
2	Insert Patient Initials or other reference associated with the patient
3	Select Patient Gender <ul style="list-style-type: none"> <li>• Male - or -</li> <li>• Female</li> </ul>
4	By clicking [Cancel] all the entered information is discarded and the application returns to the List of Patients screen (see screen N°.5 for detail)
5	Click [Save Patient] to create a new patient and move to the List of Patients screen (see screen N°.5 for detail)
	<b>Warning:</b> Under the Orthofix Terms of Use, the surgeon should never enter information that directly identifies a patient. The patient number is intended to be used as an identifying link to the patient within the surgeon's patient management system.

 **ORTHOFIX**  
Welcome Dr. TestUser


Home Page Cases Patients Account

**List of Patients** Add New Patient

Search for Patients using the boxes below.

Patient ID	Patient Initials	Gender	Date Created	Edit Patient	Delete
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<a href="#">Patient A</a>	PA	Male	27/03/2013		
<a href="#">Patient B</a>	PB	Female	27/03/2013		
<a href="#">Patient C</a>	PC	Male	27/03/2013		

Clear Search Filters


 **ORTHOFIX**  
Welcome Dr. TestUser

Home Page Cases Patients Account

**List of Cases** Add New Case

Search for Cases using the boxes below.

Patient ID	Case Number	Case Name	Anatomy	Date Created		Delete
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Patient A	001	Case 1	Left	3/27/2013	<a href="#">View</a>	
Patient A	004	Case 4	Left	3/27/2013	<a href="#">View</a>	
Patient B	002	Case 2	Right	3/27/2013	<a href="#">View</a>	
Patient C	003	Case 3	Left	3/27/2013	<a href="#">View</a>	
Patient D	004	Case 4	Left	11/6/2012	<a href="#">View or Upgrade</a>	

 **ORTHOFIX**  
Welcome Dr. TestUser

Home Page Cases Patients Account

You are attempting to open a Case that was created on an earlier software version. By clicking the Upgrade Case button, you can upgrade the case to the new version. All the data will be transferred to a new case. It is recommended that you review and confirm all of the data, because the prescription will be recalculated.

Upgrade Case Cancel




Print Page

**Case and Frame Section**



Case Viewed On: Thursday, April 04, 2013 10:00:50 AM

Case Number:	004	Patient ID:	Patient D
Case Name:	Case 4	Side Selection:	Left
Proximal Ring:	5/8 Open Posteriorly Ring - 180mm	Distal Ring:	Full Ring - 160mm
Reference Segment:	Proximal		

**5. Patient → List of Patients**

	Description
1	Click [Add New Patient] to create a new patient (see screen N°.4 for detail)
2	Select the Patient ID to lead to the list of cases associated with this patient
3	Click the icon  to edit patient details
4	Click the icon  to delete definitively a Patient. The application asks to confirm the deletion. <b>Delete action cannot be undone and implies that all the cases related to that patient will be deleted</b>
5	Use the box related to each header to filter and search the content based on a specified criteria
6	For Date Created header it is possible to set complex filter using the available function, by clicking icon 
7	To remove any filter, click [Clear search filters]

**6. Cases → List of Cases**

	Description
1	Click [Add New Case] to create a new case (see screen N°.7 for detail)
2	By default, all the cases are sorted by Patient ID. Click any of the headers (i.e. Patient ID, Case Number, Case Name, Anatomy and Date created) to change the sort
3	Click the icon  to delete definitively a Case. The application asks to confirm the deletion. <b>Delete action cannot be undone</b>
4	Use the box related to each header to filter and search the content based on a specified criteria
5	For Date Created header, it is possible to set complex filter using the available function, by clicking the 
6	The cases created with SW version 1.2 have an additional link "View" that allows a read-only version of the case to be displayed
7	In the case created with the older software version (1.1), the additional link field is VIEW OR UPGRADE


	Description
1	Click [Upgrade Case] to generate a new case using the old case input data
2	Click [Cancel] to close the page and go back to List of Cases page
3	Click [Print the page] to Print the page


The screenshot displays the TL-HEX software interface. At the top left is the TL-HEX logo with the text 'TRUELOK HEXAPOD SYSTEM'. At the top right is the ORTHOFIX logo. Below the logos is a navigation bar with 'Home Page', 'Cases', 'Patients', and 'Account'. The 'Cases' tab is active, showing a sub-menu with 'Case Data', 'Deformity Parameters', 'Frame Parameters', 'Postoperative', 'End of Correction', 'Schedule', 'Prescription', 'Report', and 'Checkup'. The 'Case Data' tab is selected, showing the following fields:

- Patient ID: Patient A (dropdown menu)
- Case Number: .001 (text input)
- Case Name: Case 1 (text input)
- Planning Created: 29/08/2012 (calendar icon)
- Side Selection: Left (radio button), Right (radio button)
- New Notes (warning icon) (text input)
- Notes History (dropdown menu)

A warning message is displayed in a box: "Warning: You are not allowed to enter or provide any information that allows, directly or indirectly, the identification of your patient (e.g. name, birth date, address, email-address, phone number etc.). Please use only an internal confidential code to identify your patient record when using this Software."

At the bottom center is a 'Next' button. At the bottom of the page are links: Privacy Policy | EULA Policy | Cookies | Contact Us | Instructions for Use.

7. Cases → Add New Case	
	Description
1	Select the Patient from the drop-down selector. If no patients have been entered, the drop-down selector will be empty. To create a Patient before beginning a case creation see screen N°4 for details
2	Assign the Case Number
3	Specify the Case Name (reference associated with this case)
4	Select the Planning Date by clicking the calendar icon
5	Select the side <ul style="list-style-type: none"><li>• Left - or -</li><li>• Right</li></ul>
6	Specify additional notes, if any
7	Review the inserted notes, if existing
8	Click [Next] to proceed with the next steps in the treatment planning process [Deformity Parameters].
	<b>Warning:</b> Under the Orthofix Terms of Use, the surgeon should never enter information that directly identifies a patient. The patient number is intended to be used as an identifying link to the patient within the surgeon's patient management system



[\[Log Out\]](#)

[Home Page](#) | [Cases](#) | [Patients](#) | [Account](#)

Case Data

Deformity Parameters

Frame Parameters

Postoperative

End of Correction

Schedule

Prescription

Report

Checkup

Case Number: 001

Case Name: Case 1

Patient ID: Patient A

Side Selection: Left

Reference Segment: ☒ Proximal ☐ Distal

AP Plane Angular Deformity (deg)

☒ Valgus ☐ Varus

ML Plane Angular Deformity (deg)

☒ Apex Anterior ☐ Apex Posterior

Rotation (deg)

☒ External ☐ Internal

AP Plane Translation (mm)

☒ Medial ☐ Lateral

ML Plane Translation (mm)

☒ Anterior ☐ Posterior

Axial Translation (mm)


☒ Short ☐ Long

Bone Length (mm)


☒ Short ☐ Long

Update Views


AP view



ML view



Axial view



Previous


Next

[Privacy Policy](#) | [EULA Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)


**8. Deformity Parameters**

	Description
1	Summary of Case Data provided in the Case Data tab
<b>Deformity Parameters</b>	
2	<p>Choose the reference segment as</p> <ul style="list-style-type: none"> <li>• Proximal - or -</li> <li>• Distal</li> </ul> <p>Refer to the Fig. 1 for the description of the moving bone segment translation depending on proximal or distal location of the reference segment.</p> <div data-bbox="421 622 1005 1001" data-label="Image"> </div> <p>Fig. 1</p> <p>In the software diagram [10], the reference segment is colored in blue while the moving segment is colored in green</p>
<b>ANGULATION</b>	
3	<p>Insert <i>coronal plane</i> angulation (AP view) in degrees. This can be either</p> <ul style="list-style-type: none"> <li>• Valgus - or -</li> <li>• Varus,</li> </ul> <p>depending whether the distal segment is bent towards or away from the midline</p>
4	<p>Insert <i>sagittal plane</i> angulation (MP view) in degrees. It is described as</p> <ul style="list-style-type: none"> <li>• Apex Anterior (procurvatum) - or -</li> <li>• Apex Posterior (recurvatum)</li> </ul>
5	<p>Insert horizontal plane angulation (Axial view) in degrees. It is described as</p> <ul style="list-style-type: none"> <li>• External rotation - or -</li> <li>• Internal rotation</li> </ul>
<b>TRANSLATION</b>	
6	<p>Insert <i>coronal plane</i> translation (AP view) in mm. It is described as</p> <ul style="list-style-type: none"> <li>• Medial - or -</li> <li>• Lateral</li> </ul>
7	<p>Insert <i>sagittal plane</i> translation (MP view) in mm. It is described as</p> <ul style="list-style-type: none"> <li>• Anterior - or -</li> <li>• Posterior</li> </ul>
8	<p>Insert <i>axial translation</i> (Axial view) in mm as</p> <ul style="list-style-type: none"> <li>• Short, used when the moving bone segment is translated (compressed) towards the reference bone segment - or -</li> <li>• Long, used when the moving bone segment is translated (distracted) away the reference bone segment</li> </ul> <p>Bone length is a clinical paramater indicating limb length discrepancy relative to the contralateral limb.</p> <ul style="list-style-type: none"> <li>• Short - or -</li> <li>• Long</li> </ul> <p>Depending on the relevant amount of bone length discrepancy in mm.</p>
9	Click [Update Views] button at any time to refresh the display according to the parameters entered
10	<p>These are three diagrams in the software.</p> <p>AP view: corresponding to the AP X-ray of the limb</p> <p>ML view: representing the ML X-ray of the limb</p> <p>Axial view: representing the view looking either up or down the limb from the reference segment</p>
11	<p>Click on [Previous] button to return to the previous tab [Case Data]</p> <p>Click on [Next] button to move to the next tab [Frame Parameters]</p>

[Log Out!](#)



Welcome Dr. TestUser



[Home Page](#) | [Cases](#) | [Patients](#) | [Account](#)

Case Data
Deformity Parameters
Frame Parameters
Postoperative
End of Correction
Schedule
Prescription
Report
Checkup

Case Number: 001  
 Case Name: Case 1

Patient ID: Patient A  
 Side Selection: Left

Print Page

**Select External Supports**

Proximal Support - Ring Type: Full Ring  
 Distal Support - Ring Type: Full Ring

Ring Size: 180mm  
 Ring Size: 180mm

**Preoperative Mounting Parameters**

**Reference Ring Parameters**

AP Translation (mm) 0 Medial Lateral

ML Translation (mm) 0 Anterior Posterior

Frame Rotation 0

**Ring Position**

Reference Ring (mm) 50 Proximal Distal

Second Ring (mm) 50


Position Relative To Deformity Apex Osteotomy/Fracture Level


Modify Horizontal Fragments Translation (mm)
 


Medial Lateral 0

Anterior Posterior 0

Update Views

**AP view**  


**ML view**  


**Axial view**  


Total	Strut 1 (mm): 185	Strut 2 (mm): 146	Strut 3 (mm): 108	Strut 4 (mm): 103	Strut 5 (mm): 125	Strut 6 (mm): 170
<b>Size</b>	Long	Medium	Short	Short	Medium	Medium
<b>Acute</b>	0	0	15	10	8	22
<b>Gradual</b>	53	3	14	14	34	1

Previous Next

[Privacy Policy](#) | [EULA Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)



**9. Frame Parameters**

	Description
1	Summary of Case Data provided in the Case Data tab
2	Click [Print Page] to print out the complete planning when it is completed

**External Supports specification**

	Ring Type available <ul style="list-style-type: none"><li>• Full Ring</li><li>• 5/8 Open Posteriorly Ring</li><li>• 5/8 Open Medially Ring</li><li>• 5/8 Open Anteriorly Ring</li></ul>	Ring Size available <ul style="list-style-type: none"><li>• 140mm</li><li>• 160mm</li><li>• 180mm</li><li>• 200mm</li><li>• 220mm</li></ul>
3	Select proximal support ring type and size from the drop-down selectors	
4	Select distal support ring type and size from the drop-down selectors	

**MOUNTING PARAMETERS SECTION** - optional. The default assumption is that both proximal and distal external supports are perpendicular to the corresponding bone segment axis and located at 50mm distance from the point of interest.

5	Adjust the position of the reference ring in the coronal plane (AP translation) in mm <ul style="list-style-type: none"> <li>• Medial - or -</li> <li>• Lateral</li> </ul>
6	Adjust the position of the reference ring in the sagittal plane (ML translation) in mm <ul style="list-style-type: none"> <li>• Anterior - or -</li> <li>• Posterior</li> </ul>
7	Adjust the reference ring position, in the axial direction in mm, described as <ul style="list-style-type: none"> <li>• Proximal - or -</li> <li>• Distal</li> </ul> Adjust the moving ring (second ring) position in mm
8	Determine the point of interest for the reference ring position <ul style="list-style-type: none"> <li>Relative to Deformity Apex - or -</li> <li>Relative to Osteotomy/Fracture Level</li> </ul> Refer to Fig. 2 for the description of the reference ring position relative to deformity apex (a) or osteotomy/fracture level (b)




Fig. 2

If the AP or ML view translation was entered into the deformity parameters section, the osteotomy/fracture level is automatically chosen as point of interest

9	Frame Rotation field with the relevant read-only value. The software will automatically change the orientation offset of the 5/8 ring
10	Click [Update Views] button, after entering the mounting parameters, will reveal the preassembled frame construct and render a set of strut lengths that will result in the required frame
11	These are the three diagrams in the software. AP view: corresponding to the AP X-ray of the limb ML view: representing the ML X-ray of the limb Axial view: representing the view looking either up or down the limb from the reference segment
12	This area shows the strut mounting parameters for all six struts. Details for each strut are described in three fields: <ul style="list-style-type: none"> <li>• [Size]: Strut size (short, medium, long)</li> <li>• [Acute]: Acute length for the strut in mm, which is read off the scale relative to the acute length orange mark</li> <li>• [Gradual]: Gradual length for the strut in mm, which is read off the scale relative to the gradual length green mark</li> </ul>
13	Click on [Previous] button to return to the previous tab [Deformity Parameters] Click on [Next] button to move to the next tab [Postoperative]

**1** **Frame Print View**

**2** **Print Frame Parameters**

**Case and Frame Section**

Case Number: 001  
Case Name: Case 1

Proximal Ring: Full Ring - 180mm  
Reference Segment: Proximal

Patient ID: Patient A  
Side Selection: Left

Distal Ring: Full Ring - 180mm

**Mounting Parameters**

Reference Ring AP Translation (mm)	Reference Ring ML Translation (mm)	Reference Ring Position (mm)	<u>100 - Proximal</u>
<u>No Translation</u>	<u>25 - Posterior</u>	Second Ring Position (mm)	<u>80</u>

**AP view**

**ML view**

**Axial view**

Total	Strut 1 (mm): 254	Strut 2 (mm): 208	Strut 3 (mm): 162	Strut 4 (mm): 149	Strut 5 (mm): 179	Strut 6 (mm): 239
<b>Size</b>	Long	Long	Long	Medium	Long	Long
<b>Acute</b>	18	49	2	33	19	3
<b>Gradual</b>	2	78	78	33	78	2

9.1 Frame Parameters - Print Page	
	Description
1	Planning summary
2	Click [Print Frame Parameters] to print out the complete planning

[\[Log Out\]](#)



Welcome Dr. TestUser



Home Page Cases Patients Account

Case Data Deformity Parameters Frame Parameters Postoperative End of Correction Schedule Prescription Report Checkup

Case Number: 001  
Case Name: Case 1

Patient ID: Patient A  
Side Selection: Left

**Reference Ring Parameters**

AP Translation (mm)  Medial Lateral

AP Angle (deg)  Medial Side Down Medial Side Up

ML Translation (mm)  Anterior Posterior

ML Angle (deg)  Anterior Side Down Anterior Side Up

Frame Rotation  External Internal

**Reference Ring Position**

Reference Ring (mm)  Proximal Distal

Position Relative To: Deformity Apex Osteotomy/Fracture Level

Modify Horizontal Fragments Translation (mm)

Medial Lateral


Anterior Posterior

Update Views

AP view



ML view



Axial view




Total	Strut 1 (mm): 234	Strut 2 (mm): 209	Strut 3 (mm): 162	Strut 4 (mm): 149	Strut 5 (mm): 179	Strut 6 (mm): 239
Size	<input type="text" value="Long"/>	<input type="text" value="Long"/>	<input type="text" value="Long"/>	<input type="text" value="Medium"/>	<input type="text" value="Long"/>	<input type="text" value="Long"/>
Acute	<input type="text" value="18"/>	<input type="text" value="49"/>	<input type="text" value="2"/>	<input type="text" value="33"/>	<input type="text" value="19"/>	<input type="text" value="3"/>
Gradual	<input type="text" value="2"/>	<input type="text" value="78"/>	<input type="text" value="78"/>	<input type="text" value="33"/>	<input type="text" value="78"/>	<input type="text" value="2"/>

Previous Next

[Privacy Policy](#) | [EULA Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)

**10. Postoperative Assessment**

	Description
<b>1</b>	Summary of Case Data provided in the Case Data tab
<b>Reference Ring Position</b>	Refer to Fig. 3 to determine the reference ring angulation, translation and rotation parameters in AP view (a), ML view (b) and axial view (c).
	 <p>Fig. 3</p>
<b>2</b>	Adjust/Enter the position of the reference ring in the coronal plane (AP translation) in mm as translation of the center of the reference ring in relation to the longitudinal axis of the reference bone segment <ul style="list-style-type: none"> <li>• Medial - or -</li> <li>• Lateral</li> </ul>
<b>3</b>	Adjust/Enter the position of the reference ring in the sagittal plane (ML translation) in mm as translation of the center of the reference ring in relation to the longitudinal axis of the reference bone segment <ul style="list-style-type: none"> <li>• Anterior - or -</li> <li>• Posterior</li> </ul>
<b>4</b>	Adjust/Enter the reference ring position, in the axial direction in mm, described as the translation of the reference ring proximal/distal along the longitudinal axis of the reference bone segment <ul style="list-style-type: none"> <li>• Proximal - or -</li> <li>• Distal</li> </ul>
<b>5</b>	Enter the reference ring angulation on the AP view in degrees as the angle between the projection of the ring and the axis of the bone segment with the medial side of the ring. <ul style="list-style-type: none"> <li>• Medial side up - or -</li> <li>• Medial side down</li> </ul>
<b>6</b>	Enter the reference ring angulation on the ML view in degrees as the angle between the projection of the ring and the axis of bone segment with the anterior side of the ring. <ul style="list-style-type: none"> <li>• Anterior side up - or -</li> <li>• Anterior side down</li> </ul>
<b>7</b>	Specify the frame rotation relative to the longitudinal axis of the reference bone segment in degrees on the axial view, described as <ul style="list-style-type: none"> <li>• External rotation - or -</li> <li>• Internal rotation</li> </ul>
<b>8</b>	Determine the point of interest for the reference ring position. See Screen 8 - Reference 8 for more details. <ul style="list-style-type: none"> <li>• Relative to Deformity Apex - or -</li> <li>• Relative to Osteotomy/Fracture Level</li> </ul>
<b>9</b>	Click [Update Views] button at any time to refresh the display according to the parameters entered. At this point, the software will generate diagrammatic models that reflect the bone deformity and frame position on the X-ray. In case of discrepancy, the surgeon should go back and check all the variables before proceeding to the next steps.
<b>10</b>	These are the three diagrams in the software. AP view: corresponding to the AP X-ray of the limb ML view: representing the ML X-ray of the limb Axial view: representing the view looking either up or down the limb from the reference segment
<b>11</b>	This area shows the strut mounting parameters for all six struts. Details for each strut are described in three fields: <ul style="list-style-type: none"> <li>• [Size]: Strut size (short, medium, long)</li> <li>• [Acute]: Acute length for the strut in mm, which is read off the scale relative to the acute length orange mark</li> <li>• [Gradual]: Gradual length for the strut in mm, which is read off the scale relative to the gradual length green mark</li> </ul> <p>Confirm or adjust the data if necessary.</p> <p>Correct any additional errors prior to proceeding with the next steps. An error is highlighted in red and the strut size and length should be corrected.</p>
<b>12</b>	Click on [Previous] button to return to the previous tab [Frame Parameters] Click on [Next] button to move to the next tab [End of Correction]

[Log Out](#)

**TL-HEX**  
TRUELOK HEXAPOD SYSTEM

Welcome Dr. TestUser

**ORTHOFIX**

Home Page Cases Patients Account

Case Data Deformity Parameters Frame Parameters Postoperative End of Correction Schedule Prescription Report Checkup

Case Number: 001  
Case Name: Case 1

Patient ID: Patient A  
Side Selection: Left

All the corrections in this section are under/over corrections relative to biomechanically aligned segment

AP Over/Under Correction (deg) 0

Valgus Varus

AP Over/Under Translation (mm) 0

Medial Lateral

ML Over/Under Correction (deg) 0

Apex Anterior Apex Posterior

ML Over/Under Translation (mm) 0

Anterior Posterior

Over/Under Rotation (deg) 0

External Internal

Bone Length (mm) 20

Shortening Lengthening

Update Views

AP view ML view Axial view

Medial Lateral Posterior Anterior Lateral Medial Posterior

Total	Strut 1 (mm): 200	Strut 2 (mm): 210	Strut 3 (mm): 252	Strut 4 (mm): 229	Strut 5 (mm): 263	Strut 6 (mm): 212
Size	Long	Long	Long	Long	Long	Long

Previous Next

[Privacy Policy](#) | [EULA Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)

11. End of Correction	
	Description
1	Summary of Case Data provided in the Case Data tab
The software assumes that, at the end of the treatment, the bone segments should be in perfect alignment. However, surgeon could override the default position as desired	
<b>ANGULATION</b>	
2	Override (Over/Under) AP correction in degrees. This can be either <ul style="list-style-type: none"> <li>• Valgus - or -</li> <li>• Varus</li> </ul>
3	Override (Over/Under) ML correction in degrees. This can be either <ul style="list-style-type: none"> <li>• Apex anterior - or -</li> <li>• Apex posterior</li> </ul>
4	Override (Over/Under) rotation in degree . It is described as <ul style="list-style-type: none"> <li>• External - or -</li> <li>• Internal</li> </ul>
<b>TRANSLATION</b>	
5	Override (Over/Under) AP translation in mm. It is described as <ul style="list-style-type: none"> <li>• Medial - or -</li> <li>• Lateral</li> </ul>
6	Override (Over/Under) ML translation in mm. It is described as <ul style="list-style-type: none"> <li>• Anterior - or -</li> <li>• Posterior</li> </ul>
7	Select the bone length <ul style="list-style-type: none"> <li>• Shortening - or -</li> <li>• Lengthening</li> </ul>
8	Click [Update Views] button at any time to refresh the display according to the parameters entered and check the bottom of the screen for struts that are out of range.
9	These are the three diagrams in the software. AP view: corresponding to the AP X-ray of the limb ML view: representing the ML X-ray of the limb Axial view: representing the view looking either up or down the limb from the reference segment
10	This area shows the strut size for all six struts.
11	Click on [Previous] button to return to the previous tab [Postoperative] Click on [Next] button to move to the next tab [Schedule]

[Log Out](#)

**TL-HEX**  
TRUELOK HEXAPOD SYSTEM

Welcome Dr. TestUser

**ORTHOFIX**

Home Page Cases Patients Account

Case Data Deformity Parameters Frame Parameters Postoperative End of Correction **Schedule** Prescription Report Checkup

Case Number: 001  
Case Name: Case 1

Patient ID: Patient A  
Side Selection: Left

Calculate

Calculate By:

☒ Daily Correction Rate (mm/day)

☐ Angular Max Speed (deg/day)

☐ Rotate Max Speed (deg/day)

☐ Days of Correction

Surgery Date

Latency Period (days)

Treatment Start Date lunedì 1 aprile 2013

Correction Time(s):

<input type="checkbox"/> 00:00:00	<input type="checkbox"/> 06:00:00	<input type="checkbox"/> 12:00:00	<input type="checkbox"/> 18:00:00
<input type="checkbox"/> 01:00:00	<input type="checkbox"/> 07:00:00	<input type="checkbox"/> 13:00:00	<input type="checkbox"/> 19:00:00
<input type="checkbox"/> 02:00:00	<input checked="" type="checkbox"/> 08:00:00	<input type="checkbox"/> 14:00:00	<input checked="" type="checkbox"/> 20:00:00
<input type="checkbox"/> 03:00:00	<input type="checkbox"/> 09:00:00	<input type="checkbox"/> 15:00:00	<input type="checkbox"/> 21:00:00
<input type="checkbox"/> 04:00:00	<input type="checkbox"/> 10:00:00	<input type="checkbox"/> 16:00:00	<input type="checkbox"/> 22:00:00
<input type="checkbox"/> 05:00:00	<input type="checkbox"/> 11:00:00	<input type="checkbox"/> 17:00:00	<input type="checkbox"/> 23:00:00



Previous Next

[Privacy Policy](#) | [EUJA Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)



12. Schedule	
	Description
1	Summary of Case Data provided in the Case Data tab
2	Specify the maximum rate of the bone segment translation (mm/day)
3	Specify the maximum rate of the bone segment rotation (degrees/day)
4	Specify the maximum rate of the bone segment angular correction (degrees/day)
5	Specified the desired days of correction
6	Indicate the surgery date
7	Indicate the latency period in days. By default this is 5 days
8	The software determines the treatment start date considering surgery date plus latency period
9	Indicate correction time/s to have the prescription calculated for one or more frame adjustments during each treatment day
10	Click on [Previous] button to return to the previous tab [End of Correction] Click on [Next] button to move to the next tab [Prescription]
11	Click on [Calculate] to calculate the speed if days of correction has been selected or the days of correction if at least one speed has been selected. (It is necessary to enter at least one value in the "Calculate By" box and "Correction Time(s)" box prior to clicking)

[Log Out](#)

[Home Page](#)
[Cases](#)
[Patients](#)
[Account](#)

**Case Data**

**Deformity Parameters**

**Frame Parameters**

**Postoperative**

**End of Correction**

**Schedule**

**Prescription**

**Report**

**Checkup**

Case Number: 001  
 Case Name: Case 1


Patient ID: Patient A  
 Side Selection: Left

Please review all information before completing and printing the prescription to ensure that it is accurate.

**Strut Adjustments in 'CLICKS'**


**Strut Reference Length (mm)**

No	Date-Time	Red Strut1	Orange Strut2	Yellow Strut3	Green Strut4	Blue Strut5	Purple Strut6	Red Strut1	Orange Strut2	Yellow Strut3	Green Strut4	Blue Strut5	Purple Strut6	View Image	See In Report
0	04/09/2012 00:00	0	0	0	0	0	0	14	40	35	35	35	33	<a href="#">View</a>	<a href="#">In Report</a>
1	04/09/2012 08:00	0	0	+1	+2	+1	0	14	40	34	34	34	33	<a href="#">View</a>	<a href="#">In Report</a>
2	04/09/2012 20:00	-1	+1	+2	+3	+2	0	14	40	34	34	32	33	<a href="#">View</a>	<a href="#">In Report</a>
3	05/09/2012 08:00	0	0	+1	+3	+1	0	15	39	33	31	33	33	<a href="#">View</a>	<a href="#">In Report</a>
4	05/09/2012 20:00	-1	+1	+2	+3	+2	0	15	39	32	29	32	33	<a href="#">View</a>	<a href="#">In Report</a>
5	06/09/2012 08:00	0	0	+1	+4	+1	0	15	39	31	28	31	33	<a href="#">View</a>	<a href="#">In Report</a>
6	06/09/2012 20:00	-1	+1	+2	+3	+2	0	15	39	31	26	31	33	<a href="#">View</a>	<a href="#">In Report</a>
7	07/09/2012 08:00	0	0	+1	+3	+1	0	16	38	30	25	30	33	<a href="#">View</a>	<a href="#">In Report</a>
8	07/09/2012 20:00	-1	+1	+2	+3	+2	0	16	38	29	23	29	33	<a href="#">View</a>	<a href="#">In Report</a>
9	08/09/2012 08:00	0	0	+1	+3	+1	0	16	38	28	21	28	33	<a href="#">View</a>	<a href="#">In Report</a>
10	08/09/2012 20:00	-1	+1	+2	+3	+2	0	16	38	28	20	28	33	<a href="#">View</a>	<a href="#">In Report</a>
11	09/09/2012 08:00	-1	0	+1	+3	+1	0	17	37	27	18	27	33	<a href="#">View</a>	<a href="#">In Report</a>
12	09/09/2012 20:00	0	+1	+2	+3	+2	0	17	37	26	17	26	33	<a href="#">View</a>	<a href="#">In Report</a>
13	10/09/2012 08:00	-1	0	+2	+3	+1	0	17	37	25	15	25	33	<a href="#">View</a>	<a href="#">In Report</a>




AP view

AP Angle: 25 deg, Valgus  
AP Translation: 0 mm



ML view

ML Angle: 25 deg, ApexAnterior  
ML Translation: 0 mm



Axial view

Rotation: 0 deg  
Axial Translation: 0 mm

Previous

Next


[Privacy Policy](#) | [EULA Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)

**13. Prescription**

	Description
1	Summary of Case Data provided in the Case Data tab
2	Click [Print Prescription] to generate a prescription in .pdf format. This document can be saved for the record and printed as hard copy to be issued to the patient
3	Represent the adjustment for each strut by number of clicks (1/2 rotation of the strut adjustment knob). It can be <ul style="list-style-type: none"> <li>• positive (if strut length increases) - or -</li> <li>negative (if strut length decreases)</li> </ul>
4	Represent the gradual adjustment scale value in millimeters as reference for each strut
5	Click [View] in the prescription table to have a look to the three views of the deformity and the frame for the corresponding day/time of the prescription
6	Click [In Report] to see the same adjustment row in the Report tab
7	These are the three diagrams in the software considering the deformity and the frame for the corresponding day/time of the prescription. <p>AP view: corresponding to the AP X-ray of the limb  ML view: representing the ML X-ray of the limb  Axial view: representing the view looking either up or down the limb from the reference segment</p>
8	Click on [Previous] button to return to the previous tab [Schedule] Click on [Next] button to move to the next tab [Report]

The prescription row will be highlighted when strut readjustments (shaded blue) or exchanges (shaded red) are required.


The rows of lighter shading indicate the allowable range of days that is suitable for the readjustment/exchange; the heavier shaded row indicates the last possible day for the strut change



**TL-HEX**  
TRUELOK HEXAPOD SYSTEM

Welcome Dr. TestUser

[\[Log Out\]](#)



Home Page
Cases
Patients
Account

Case Data
Deformity Parameters
Frame Parameters
Postoperative
End of Correction
Schedule
Prescription
Report
Checkup

Case Number: 001

Case Name: Case 1

Patient ID: Patient A


Side Selection: Left

[Print Report](#)


Strut Length A-Acute / G-Gradual

No	Date-Time	Red		Orange		Yellow		Green		Blue		Purple		Images	View In							
		1:	A	G	2:	A	G	3:	A	G	4:	A	G	5:	A	G	6:	A	G			
1	01/04/2013 08:00	long	0	14	long	0	40	long	0	34	med	16	34	long	11	34	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
2	01/04/2013 20:00	long	0	14	long	0	40	long	0	33	med	16	30	long	11	33	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
3	02/04/2013 08:00	long	0	15	long	0	39	long	0	32	med	16	28	long	11	32	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
4	02/04/2013 20:00	long	0	16	long	0	38	long	0	30	med	16	26	long	11	30	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
5	03/04/2013 08:00	long	0	16	long	0	38	long	0	29	med	16	23	long	11	29	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
6	03/04/2013 20:00	long	0	16	long	0	38	long	0	28	med	16	20	long	11	28	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
7	04/04/2013 08:00	long	0	17	long	0	38	long	0	26	med	16	18	long	11	26	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
8	04/04/2013 20:00	long	0	18	long	0	37	long	0	25	med	16	15	long	11	25	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
9	05/04/2013 08:00	long	0	18	long	0	36	long	0	24	med	16	12	long	11	24	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
10	05/04/2013 20:00	long	0	18	long	0	36	long	0	22	med	16	10	long	11	22	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
11	06/04/2013 08:00	long	0	19	long	0	36	long	0	22	med	16	8	long	11	22	long	15	33	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
12	06/04/2013 20:00	long	0	20	long	0	36	long	0	20	med	16	5	long	11	20	long	15	34	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
13	07/04/2013 08:00	long	0	20	long	0	35	long	0	18	med	16	3	long	11	19	long	15	34	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
14	07/04/2013 20:00	long	0	20	long	0	35	long	0	18	med	16	0	long	11	18	long	15	34	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
15	08/04/2013 08:00	long	0	21	long	0	34	long	0	16	long	6	77	long	11	16	long	15	34	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
16	08/04/2013 20:00	long	0	22	long	0	34	long	0	15	long	6	74	long	11	15	long	15	34	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>
17	09/04/2013 08:00	long	0	22	long	0	34	long	0	14	long	6	72	long	11	14	long	15	34	<a href="#">View</a>	<a href="#">Prescription</a>	<a href="#">Checkup</a>


AP view



ML view



Axial view



AP Angle: 25 deg, Valgus

AP Translation: 0 mm

ML Angle: 25 deg, ApexAnterior

ML Translation: 0 mm

Rotation: 0 deg

Axial Translation: 0 mm


Total	Strut 1 (mm): 224	Strut 2 (mm): 198	Strut 3 (mm): 203	Strut 4 (mm): 130	Strut 5 (mm): 214	Strut 6 (mm): 220
<b>Size</b>	Long	Long	Long	Medium	Long	Long
<b>Acute</b>	0	0	0	16	11	15
<b>Gradual</b>	14	40	35	35	35	33

[Previous](#) [Next](#)


[Privacy Policy](#) | [EULA Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)

14. Report	
	Description
1	Summary of Case Data provided in the Case Data tab
2	Click [Print Report] generate the report in .pdf format. This document can be saved for the record or printed as hard copy
3	Provides a more detailed prescription. In addition to the patient prescription it includes <ul style="list-style-type: none"> <li>• [St] Strut size</li> <li>• [A] Acute adjustment in millimeters</li> <li>• [G] Gradual adjustment in millimeters for each strut</li> </ul>
4	Click [View] in the report table to have a look to the three views of the deformity and the frame for the corresponding day/time of the prescription and the struts parameters and length on the desired day
5	Click [Prescription] in the report table to see the same adjustment row in the Prescription tab
6	Click [Checkup] in the Report table to launch the Checkup for the selected day and time
7	These are the three diagrams in the software considering the deformity and the frame for the corresponding day/time of the prescription. <p>AP view: corresponding to the AP X-ray of the limb  ML view: representing the ML X-ray of the limb  Axial view: representing the view looking either up or down the limb from the reference segment</p>
8	Click on [Previous] button to return to the previous tab [Prescription] Click on [Next] button to move to the next tab [Checkup]
<p>The prescription row will be highlighted when strut readjustments (shaded blue) or exchanges (shaded red) are required.  The rows of lighter shading indicate the allowable range of days that is suitable for the readjustment/exchange; the heavier shaded row indicates the last possible day for the struct change</p>	

[Log Out](#)



Welcome Dr. TestUser



[Home Page](#) | [Cases](#) | [Patients](#) | [Account](#)

Case Data
Deformity Parameters
Frame Parameters
Postoperative
End of Correction
Schedule
Prescription
Report
Checkup

Case Number: 001

Case Name: Case 1


Patient ID: Patient A

Side Selection: Left

Treatment Date:


Displaying prescription line No. 0 based on treatment range of : 01/04/2013 - 21/04/2013

AP view




Medial      Lateral

ML view



Posterior      Anterior

Axial view



Lateral      Medial

Posterior

AP Angle: 25 deg, Valgus

AP Translation: 0 mm

ML Angle: 25 deg, ApexAnterior

ML Translation: 0 mm

Rotation: 0 deg

Axial Translation: 0 mm

Total	Strut 1 (mm): 224	Strut 2 (mm): 198	Strut 3 (mm): 203	Strut 4 (mm): 130	Strut 5 (mm): 214	Strut 6 (mm): 220
<b>Size</b>	Long	Long	Long	Medium	Long	Long
<b>Acute</b>	0	0	0	16	11	15
<b>Gradual</b>	14	40	35	35	35	33

To Create New Case

New Case Number:

New Case Name:

Create New Case

[Previous](#)

[Privacy Policy](#) | [F.I.B.A. Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)

1

2

6

3

4

5

7

15. Checkup	
	Description
1	Summary of Case Data provided in the Case Data tab
	Provides the position of the bone segments and the frame with corresponding strut adjustment values at any particular day of treatment (deformity correction)
2	Select the treatment date. As default, the screen shows the data for the current day
3	Enter new case detail - Case Number
4	Enter new case detail - Case Name
5	Click [Create New Case] to open a new Case data screen for the newly generated case
6	These are the three diagrams in the software considering the deformity and the frame for the corresponding day/time of the prescription.  AP view: corresponding to the AP X-ray of the limb ML view: representing the ML X-ray of the limb Axial view: representing the view looking either up or down the limb from the reference segment  All the deformity and frame parameters are transferred from the previous case at the date of check-up
7	Click on [Previous] button to return to the previous tab [Report]
	Select the treatment date. As a default, the screen shows the data for the current day. The patient, based on the starting point as chosen from the checkup screen

**TL-HEX**  
TRUELOK HEXAPOD SYSTEM

[Log Out](#)

**ORTHOFIX**

Home Page Cases Patients **Account**  
Change Password

**Change Password**

Account Information

Old Password:

New Password:

Confirm New Password:

Cancel Change Password

[Privacy Policy](#) | [EULA Policy](#) | [Cookies](#) | [Contact Us](#) | [Instructions for Use](#)

1

2

3

4

5



16. Change Password	
	Description
1	Type the old password
2	Chose and type the new password. It must be 6 or more characters
3	Confirm by typing the new password again. The software will validate that the two passwords match
4	Click [Cancel] to cancel the action and come back to the home page
5	Click [Change Password] to submit the password change

Manufactured by:  
ORTHOFIX Srl  
Via Delle Nazioni 9  
37012 Bussolengo (Verona)  
Italy

Telephone +39 045 6719000  
Fax +39 045 6719380



Your Distributor is: