| | | | |
|----|--|----|--|
| EN | TrueLok WIRE TENSIONER (54-1139) 2 | ZH | TrueLok 张力钳 (54-1139) 28 |
| IT | PINZA TENDIFILO TrueLok (54-1139) 4 | CS | NAPÍNAČ DRÁTU TrueLok (54-1139) 30 |
| FR | TENDEUR POUR BROCHE TrueLok (54-1139) 6 | PL | NAPREŽACZ DRUTÓW TrueLok (54-1139) 32 |
| DE | TrueLok DRAHTSPANNER (54-1139) 8 | SL | NAPENJALNIK ŽICE TrueLok (54-1139) 34 |
| ES | TENSOR DE ALAMBRE TrueLok (54-1139) 10 | KO | TrueLok 와이어 텐션너 (54-1139) 36 |
| DA | TrueLok-TRÅDSTRAMMER (54-1139) 12 | AR | آداة شد السلك (54-1139) 38 |
| FI | TrueLok-PIIKINKIRISTIN (54-1139) 14 | BG | Обтегач за проводници TrueLok WIRE TENSIONER (54-1139) 40 |
| NO | TrueLok VAIERSTRAMMER (54-1139) 16 | LT | „TrueLok“ VIELOS ĮTEMPIKLIS (54-1139) 42 |
| NL | TrueLok-SPANAPPARAAT (54-1139) 18 | MS | PENEGANG WAYAR TrueLok (54-1139) 44 |
| PT | TENSIONADOR DE FIO TrueLok (54-1139) 20 | RO | DISPOZITIV DE TENSIONARE A FIRELOR TrueLok (54-1139) 46 |
| SV | TrueLok TRÅDSPÄNNARE (54-1139) 22 | TR | TrueLok KABLO GERİCİ (54-1139) 48 |
| EL | ENTATHRAS ΣΥΡΜΑΤΟΣ TrueLok (54-1139) 24 | HR | ZATEZAČ ŽICE TrueLok (54-1139) 50 |
| JA | TrueLok ワイヤーテンショナー (54-1139) 26 | LV | TrueLok STIEPĻU SPRIEGOTĀJS (54-1139) 52 |

ORTHOFIX Srl
 Via Delle Nazioni 9 - 37012 Bussolengo (Verona) - Italy
 Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

0123

PQWTN F 03/21 (0421217)

Instructions For Use are subject to change; the most current version of each Instruction For Use is always available online.

Important information - please read prior to use. See also instruction leaflet PQTLK for the TrueLok™ Systems and PQRMD for reusable medical devices.

Keep this leaflet: you may need to read it again

TrueLok WIRE TENSIONER (54-1139)

Manufacturer name

ORTHOFIX SRL, via Delle Nazioni 9, 37012 Bussolengo (VR) Italy
Tel. 0039 (0) 45 6719000 - Fax 0039 (0) 45 6719380

- 1) The Wire Tensioner (54-1139) allows tensioning wires up to 130kg. An indication of the applied tension can be read from a scale (fig. 1).


Figure 1

- 2) Only when tensioning 1.5mm diameter wires, apply the specific Tension Limiter (54-1141) to physically limit the applied tension to 90kg.

Make sure to insert the Tension Limiter correctly, paying attention to the laser mark (fig. 2).


Figure 2

- 3) To exchange the different types of Wire Tensioner Tips, remove the Wire Tensioner's Tip by loosening the locking set-screw with the hex-driver (52-1018) (fig. 3). Do not remove the locking set-screw for any reason.


Figure 3

| | | | |
|----------|---------|---------|---|
| 20116735 | 54-1151 | 54-1152 | ✓ |
| 20116736 | 54-1251 | 54-1252 | ✓ |
| | 54-1151 | 54-1152 | ✓ |
| 20116731 | 54-1251 | 54-1252 | ✗ |

- 4) Using appropriate tools thoroughly clean the whole cannulation (fig. 4) to remove any debris that might be present. Additional information on cleaning and reprocessing can be found in the leaflet PQRMD.


Figure 4

- 5) Apply one drop of medical grade lubricant on moving elements when needed.

NOTE: Products that show excessive fading of marked product code and lot, thus preventing clear identification and traceability, must NOT BE USED.

All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

| Symbol | Description | |
|---|--|---|
| MD | Medical Device | |
|   | Consult instructions for use or consult electronic instructions for use | Caution: Consult instructions for use for important cautionary information |
|  | Single Use. Do not re-use | Orthofix note: discard appropriately after the use (treatment) on patient |
| STERILE R | Sterile. Sterilised by irradiation | |
|  | Non-sterile | |
|  | Double sterile barrier system | |
| REF | LOT | Catalogue number |
|  | Use-by date (year-month-day) | |
| CE | CE 0123 | CE marking in conformity to applicable European Medical Device Directives/Regulations |
|  | Date of manufacture | Manufacture |
|  | Do not use if package is damaged and consult instructions for use | |
| Rx Only | Caution: Federal law (USA) restricts this device to sale by or on the order of a physician | |

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