

# L. Contours Proximal Humeral Plate (PHP) + LSP Titanium Humeral Plate



1

## 1 Contours PHP

See Operative Technique PH-1201-OPT  
"Contours PHP Proximal Humeral Plate".

CE 0123



Manufactured by:  
**ORTHOFIX Srl**  
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy  
Telephone +39 045 6719000, Fax +39 045 6719380

# 1. Contours PHP Components



**1**  
**Contours PHP (Sterile)**  
 99-TP3001 Left  
 99-TP3000 Right


**2**  
**Diaphyseal Screw  $\varnothing$  4.5mm (Sterile)**  
 99-TP1422 Length 22mm  
 99-TP1424 Length 24mm  
 99-TP1426 Length 26mm  
 99-TP1428 Length 28mm  
 99-TP1430 Length 30mm  
 99-TP1432 Length 32mm  
 99-TP1434 Length 34mm  
 99-TP1436 Length 36mm

**3**  
**Revision Screw  $\varnothing$  6.5mm (Sterile)**  
 99-TP1526 Length 26mm  
 99-TP1528 Length 28mm  
 99-TP1530 Length 30mm  
 99-TP1532 Length 32mm  
 99-TP1534 Length 34mm

**4**


**Main Locking Screw  $\varnothing$  6.0mm (Sterile)**  
 99-TP1640 Length 40mm  
 99-TP1645 Length 45mm  
 99-TP1650 Length 50mm  
 99-TP1655 Length 55mm  
 99-TP1660 Length 60mm

**5**  
**Fine Threaded Screw  $\varnothing$  2.2mm, pack of 2 (Sterile)**  
 99-TP1720 Length 120/20mm  
 99-TP1725 Length 120/25mm  
 99-TP1730 Length 120/30mm  
 99-TP1735 Length 120/35mm  
 99-TP1740 Length 120/40mm  
 99-TP1745 Length 120/45mm  
 99-TP1750 Length 120/50mm  
 99-TP1755 Length 120/55mm



## 2. Contours PHP Instrumentation (Sterile)

(Out of Tray)




**1**  
**99-TP1111** K-Wire (Ø 2mm) L. 200mm  
(pack of 4)

**2**  
**99-1-TP3107** Cannulated Drill Bit (Ø 4.8mm)  
quick connect (sterile)

**3**  
**99-1-TP1104** Drill Bit (Ø 3.2mm) quick  
connect (sterile)



### 3. Sterilization Box Contours PHP



1

#### 1 TP3995 Contours PHP Sterilization Box, empty

Can accommodate:

- 1xTP3106 Quick connect Screwdriver for hexagonal head
- 1xTP1145 Ruler
- 1xTP3102 3.5mm quick connect Screwdriver
- 1x1-TP3107 Cannulated Drill Bit (Ø 4.8mm) quick connect
- 2x1-TP1104 Drill Bit (Ø 3.2mm) quick connect
- 1xTP3150 Handle (soft touch) quick connect
- 1xTP3108 Plate Holder
- 1xTP1100 Drill Guide (Ø 3.2mm)
- 1xTP1101 Locking T-Wrench

