


O. Radiolucent Wrist Fixator



1

1 99-36501 RadioLucent Wrist Fixator Kit, sterile

Consisting of: 1 RadioLucent Wrist Fixator, 1 Compression-Distraction Unit, 1 Allen Wrench 4mm (36017), 2 Drill Guides, 1 Drill Bit \varnothing 2.7mm (1355001), 1 Screw Template, 1 T-Wrench (M210), 2 Self-drilling Screws, shaft \varnothing 4mm, thread \varnothing 3.3-3mm, total length 70mm, thread length 20mm, 2 Self-drilling Screws, shaft \varnothing 4mm, thread \varnothing 3.3-3mm, total length 80mm, thread length 25mm and 4 Bone Screw Covers.

See Quick Reference Guide PG 150 "The RadioLucent Wrist Fixator" for application instructions.



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380