



QUICK REFERENCE GUIDE


See Operative Technique HF-1501-OPT "Orthofix Chimaera Hip Fracture System Trochanteric Nailing System"




1 Opening of the medullary canal




2 Nail insertion




3 Lateral cortex predrilling




4 3,2mm guide wire insertion




5 Trocar removal



6 Main lag screw measurements




7 Lag screw gradual reamer setting




8 Main lag screw reaming




9 Main lag screw connection




10 Main lag screw insertion




11 OPTIONAL STEP
Compression of the sliding lag screw




12 OPTIONAL STEP
Insertion of the supplementary lag screw



13 Distal locking screw sizing



14 Distal locking screw insertion



15 End cap insertion