

Chimaera

Hip Fracture System™

Innovation
takes root
in great
standards

Chimaera Hip
Fracture System™
designed to
leverage decades
of **knowledge**
and **experience** in
hip nailing

 **ORTHOFIX**®
ORTHOPEDICS

Chimaera

Hip Fracture System™

Indicated for treatment of **stable and unstable** pertrochanteric, intertrochanteric and subtrochanteric **fractures**

All-in-one instrument tray
Intuitive **color coded**
instrument selection

Instrumentation allows
reproducible, easy and
accurate measurement and may
reduce related x-ray exposure

Supplementary lag screw
anti-rotational stability

Self-locking
sliding lag screw

Designed to
reduce lateral pain during
weight-bearing and the
related **reoperation rate**

TECHNICAL DETAILS	
Short nails length:	180mm
Long nails length:	280-460mm (20mm increments)
Angulations:	125°, 130°
Metaphyseal angle:	5°
Proximal diameter:	15.5mm
Distal diameter:	10mm, 11mm
Cannulated guided end caps	



1:1 scale



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380

CE 0123

orthofix.com

HF-1705-PL-E0 A 10/22

Distributed by:



 **ORTHOFIX®**
ORTHOPEDICS